

# NIPD-RhD – Rhesus test

## Fact sheet

- Non-invasive prenatal RhD genotyping from maternal blood
- For RhD negative pregnant women with singleton pregnancy
- Between gestational weeks 11 + 0 and 25 + 0 (recommended)
- CE-marked IVD<sup>1,2</sup> – reimbursed in France
- Test performed by Eurofins Biomnis, France
- Delivery time 10 working days per analysis (Monday to Friday, without holidays)

**100%<sup>3</sup>**

**Sensitivity**

**95.1%<sup>3</sup>**

**Specificity**

### Limitations

- False-positive results are possible in case of a rare genotype, non-functional RhD variant<sup>4</sup> or bone marrow donation<sup>5</sup>
- A test result cannot be reported if the mother is carrier of the RhD gene
- In general, false-negative results are possible<sup>6,7</sup>

### Exemplary diagnostic strategy for RhD negative pregnant women

#### Non-immunised patient

**Blood sampling & express delivery**  
(responsible doctor)

**Test result**  
(generally 10 working days)


#### Immunised patient

**Blood sampling & express delivery**  
(responsible doctor)

**Test result**  
(generally 10 working days)

### Timeframe for blood sampling

We recommend blood sampling latest in gestational week 25 + 0. This will ensure a test result latest in gestational week 29 – including a repetition, if necessary.



### Repeat analysis to confirm a negative test result

- You will receive a special blood collection set for the repeat analysis.
- Fill both blood tubes
- Tick „repeat blood sample (if 1st analysis failed)“ on the form
- Send both blood samples together with the completed form

\* due to biological limitations

1 Rouillac-Le Sciellour C et al. (2007). Noninvasive fetal RHD genotyping from maternal plasma. Use of a new developed Free DNA Fetal Kit RhDR. *Transfus Clin Biol*, 2007 Dec;14(6):572-7. doi: 10.1016/j.tracl.2008.01.003.

2 Mackie FL et al. (2017). The accuracy of cell-free fetal DNA-based non-invasive prenatal testing in singleton pregnancies: a systematic review and bivariate meta-analysis. *BJOG*, 2017 Jan;124(1): 32-46. doi: 10.1111/1471-0528.14050.

3 Instruction manual: Institut de Biotechnologies Jaques Boy. Free DNA Fetal KitR RhD. Noninvasive fetal RHD genotyping from free fetal DNA in maternal RhD-Negative pregnant women blood (Real-Time PCR). Product identification: 502080233. Manual version 16/03/2018.

4 Flegel WA (2007). Genetik des Rhesus-Blutgruppensystems. *Dtsch Arztebl* 2007, 104(10): A-651-657/B-573/C-549.

5 Thurik FF et al. (2016). Fetal RHD genotyping after bone marrow transplantation. *Transfusion*, 56: 2122-2126. doi:10.1111/trf.13669.

6 Clausen FB et al. (2014). Routine noninvasive prenatal screening for fetal RHD in plasma of RhD-negative pregnant women – 2 years of screening experience from Denmark. *Prenat Diagn*, 34:1000-1005. doi: 10.1002/pd.4419.

7 de Haas, M et al. (2012). A nation-wide fetal RHD screening programme for targeted antenatal and postnatal anti-D. *ISBT Science Series*, 7: 164-167. doi:10.1111/j.1751-2824.2012.01600.x.

# 22q11.2 microdeletion

## Fact sheet

- Microdeletion (loss of genetic information) of about 3 million bases (3 megabases) on the long arm of chromosome 22 in position 11
- Associated with DiGeorge syndrome and velo-cardio-facial syndrome
- Occurs spontaneously in more than 90% of those affected, can also be inherited in 6% – 28%<sup>1</sup> of the cases
- Clinical abnormalities, particularly heart defects, may have different manifestations, depending on the severity of the disease

### Determination of a 22q11.2 microdeletion

**Non-invasive prenatal test to exclude or detect a 22q11.2 microdeletion<sup>2</sup>** for singleton pregnancy (using the proprietary & CE-marked *dmap* software)

**NEW** – can be selected in combination with both PrenaTest<sup>®</sup> options 2 and 3

Test result provided only with **cffDNA level  $\geq$  11%** in addition to other quality criteria

If quality criteria are not met: no repetition/no invoicing

**NEW** – post analysis possible with PrenaTest<sup>®</sup> data  $\leq$  3 months; without new blood sample

**Test results require 2 additional working days<sup>3</sup>** after analysis of test options 2/3 – to be invoiced separately

### Test accuracy of the PrenaTest<sup>®</sup> for the determination of the 22q11.2 microdeletion

**Phase 1** – A total of 469 samples were tested, of which 175 (37.3%) met the quality criteria. Three positive samples with a 22q11.2 microdeletion were correctly determined (3/3, 100%). All negative samples were correctly classified (172/172, 100%). There were no false-positive or false-negative results.

**Phase 2** – In a final internal blinded study, 20 samples from Phase 1 were examined. All samples were classified correctly. Due to the low number of cases a concrete test sensitivity and specificity cannot be derived.

### Limitations

#### cffDNA level

The detection or exclusion of a 22q11.2 microdeletion is currently only possible for specimens which have a cffDNA level of at least 11% (in addition to other quality criteria).

#### Mother is the carrier of the microdeletion

It is possible that the mother is the carrier of the 22q11.2 microdeletion. This can lead to discordant (false-positive) test results.

#### Size of the microdeletion

In more than 85% of affected persons, the deletion includes a region measuring approx. 2.5 megabases in the 22q11.2 region of chromosome 22. This region is investigated with the PrenaTest<sup>®</sup>. A small percentage of affected persons has an even smaller deletion or point mutation in the affected region which cannot be detected with the PrenaTest<sup>®</sup>. This can lead to discordant (false-negative) test results.

#### Microduplication in the investigated gene region

In the investigated gene region of chromosome 22, microduplications can also occur. These are not determined with the PrenaTest<sup>®</sup>.

### Reasons for the genetic examination

This examination is particularly useful in cases of abnormalities on ultrasound, for example during organ screening, which may correlate with a DiGeorge or velo-cardio-facial syndrome:

- Congenital heart defect
- Detection of an aberrant subclavian artery
- Kidney malformation
- Growth retardation
- Cleft lip and palate
- Increased nuchal translucency
- Small thymus

### Clinical lab routine

- Over 1,000 tests successfully performed since June 2016
- So far no discordant results based on feedback from doctors
- But: false-negative or false-positive results cannot be excluded

### Please note:

The cffDNA level is highly individual and increases in later gestational weeks. Therefore it is possible that due to a cffDNA level that is too low, no test result can be achieved. A new blood sample is not requested provided the cffDNA level is sufficient for determining fetal trisomies 13, 18, 21 as well as gonosomal aneuploidies.

Week 12 – 14	Week 19 – 23
Samples with cffDNA $\geq$ 11%	Samples with cffDNA $\geq$ 11%
47.25%	52.78%
Samples with cffDNA $<$ 11%	Samples with cffDNA $<$ 11%
53.75%	47.22%

Correlation between cffDNA level and gestational week based on data of lab routine (12/2018 – 03/2019)

1 Fernández L et al. Higher frequency of uncommon 1.5-2 Mb deletions found in familial cases of 22q11.2 deletion syndrome.

Am J Med Genet A. 2005 Jul 1;136(1):71-5.

2 OMIM# of the disease: 188400 (DGS), 192430 (VCFS)

3 Monday to Friday, except holidays