

Methods and diagnostic quality

In the NIPD-RhD, cell-free DNA from the placenta is tested. Using qPCR, the exons 5, 7 and 10 of the RhD gene are amplified and evaluated. If the test result is positive for all three exons, the analyzed sample is classified as RhD-positive. If the test result is positive for only one or two exons, the overall result is evaluated as being not diagnostically conclusive. All test results of the three exons must be negative in order for the overall result of the sample to be classified as RhD-negative.

Test result			Classification
Exon 5	Exon 7	Exon 10	
+	+	+	RhD-positive
-	-	-	RhD-negative
1 or 2 exons positive			not diagnostically conclusive

IQWiG: Prenatal test is equivalent to the postnatal test

In its final report (March 2018)¹⁰, the Institute for Quality and Efficiency in Health Care (IQWiG) certifies that the non-invasive prenatal determination of the fetal Rhesus factor has a very high level of sensitivity and specificity. Because of the high level of sensitivity of the test, only very few cases of anti-D prophylaxis which were erroneously not administered antepartum can be expected. In addition, the postnatal tests also yielded false-negative results at rates similar to those of the prenatal tests. The prenatal test is therefore equivalent to the postnatal test.

- 1 Flegel WA (2007). Genetik des Rhesus-Blutgruppensystems [Genetics of the Rhesus blood group system]. Dtsch Arztebl 2007, 104(10): A-651-657/B-573/C-549.
- 2 Hemotherapy guideline, overall amendment 2017, chapter 4.12.1.5
- 3 237th statement from the German Society of Gynecology and Obstetrics (DGGG). Nichtinvasive Bestimmung des fetalen Rhesusfaktors zur Vermeidung einer mütterlichen Rhesus-Sensibilisierung im Rahmen der Vorsorgeuntersuchung gemäß Mutterschafts-Richtlinien (Mu-RL) [Non-invasive determination of the fetal Rhesus factor to avoid maternal Rhesus sensitization during screening in accordance with the maternity directive (Mu-RL)]. www.dggg.de/leitlinien-stellungnahmen/stellungnahmen/nichtinvasive-bestimmung-des-fetalen-rhesusfaktors-413, accessed on 01/21/2019.
- 4 Thuri FF et al. (2016). Fetal RHD genotyping after bone marrow transplantation. Transfusion, 56: 2122-2126. doi:10.1111/trf.13669.
- 5 Clausen FB et al. (2014). Routine noninvasive prenatal screening for fetal RHD in plasma of RhD-negative pregnant women – 2 years of screening experience from Denmark. Prenat Diagn, 34:1000-1005. doi: 10.1002/pd.4419.
- 6 de Haas, M et al. (2012). A nation wide fetal RHD screening programme for targeted antenatal and postnatal anti D. ISBT Science Series, 7: 164-167. doi:10.1111/j.1751-2824.2012.01600.x.
- 7 Instruction manual: Jaques Boy Biotechnology Institute. Free DNA Fetal KitR RhD. Noninvasive fetal RHD genotyping from free fetal DNA in maternal RhD-Negative pregnant women blood (Real-Time PCR). Product identification: 502080233. Manual version 16/03/2018.
- 8 Rouillac-Le Sciellour C et al. (2007). Noninvasive fetal RHD genotyping from maternal plasma. Use of a new developed Free DNA Fetal Kit RhDR. Transfus Clin Biol, 2007 Dec;14(6):572-7. doi: 10.1016/j.traccl.2008.01.003.
- 9 Mackie FL et al. (2017). The accuracy of cell-free fetal DNA-based non-invasive prenatal testing in singleton pregnancies: a systematic review and bivariate meta-analysis. BJOG, 2017 Jan;124(1): 32-46. doi: 10.1111/1471-0528.14050.
- 10 IQWiG report no. 607. Nichtinvasive Bestimmung des fetalen Rhesusfaktors zur Vermeidung einer mütterlichen Rhesussensibilisierung [Non-invasive determination of the fetal Rhesus factor to avoid maternal Rhesus sensitization]. Order D16-01. Final report, version 1.0., www.iqwig.de/download/D16-01_Bestimmung_fetalen-Rhesusfaktor_Abschlussbericht_V1-0.pdf. Status: 03/20/2018, accessed on 08/14/2018.

Eurofins LifeCodexx First NIPT provider in Europe

We have been developing clinically validated, non-invasive prenatal tests (NIPT) since 2010. In 2012, the PrenaTest® was launched on the market as Europe's first NIPT. It has now become firmly established in many gynecological practices in Europe, the Middle East, and Asia as a reliable, fast and safe test method. Since the start of 2018, LifeCodexx has belonged to Eurofins Scientific, a leading international laboratory group, and has strengthened the technology portfolio of the Clinical Diagnostics business unit.

Expert hotline: +49 (0) 7531-97694844

www.lifecodexx.com

Eurofins LifeCodexx GmbH

Line-Eid-Straße 3, 78467 Konstanz, Germany

Phone +49 (0) 7531-9769460, Fax +49 (0) 7531-9769480, info@lifecodexx.com

© LifeCodexx / PraenaTest / PrenaTest are registered trademarks of Eurofins LifeCodexx GmbH, Germany



Eurofins LifeCodexx is part of



First-class medical care, available worldwide:

- Over 65,000 patients per day
- Over 110 million tests per year
- Over 800 sample collection centers
- Over 250 city laboratories
- 60 on-site hospital laboratories
- From sample collection to medical counseling
- Clinical interpretation provided
- Accredited laboratories

www.eurofins.com

Follow us:

- @EurofinsGroup
- www.facebook.com/EurofinsGroup
- www.linkedin.com/company/Eurofins



NIPD-RhD

Non-invasive prenatal
RhD genotyping

lifecodexx Information for physicians

03



NIPD-RhD – The Rhesus test from Eurofins LifeCodexx

Anti-D prophylaxis often unnecessary

In general, anti-D prophylaxis is recommended to RhD-negative pregnant women in the 28th week of pregnancy, regardless of whether their unborn child or the child's father is RhD-positive. However, the prophylaxis is unnecessary in about 40% of the RhD-negative women since their unborn children are also RhD-negative.¹

Non-invasive prenatal RhD determination from maternal blood

Now you have the option of determining the fetal Rhesus factor non-invasively from the blood of the RhD-negative mother (NIPD-RhD), with a high level of diagnostic reliability. As a result, you can decide whether anti-D prophylaxis is indicated in order to use it specifically only for those RhD-negative pregnant women who are actually expecting an RhD-positive child.

Anti-D prophylaxis in pregnant women is not necessary if the fetus was determined to be RhD-negative using a validated method.

After delivery, the RhD characteristic is to be determined, preferably from umbilical cord blood.

Hemotherapy guideline, overall amendment 2017²

Benefits and use

Benefits of the NIPD-RhD

According to the German Society of Gynecology and Obstetrics (DGGG), you avoid the following adverse effects in those pregnant women in whom prophylaxis is not indicated³:

- Potential infection through the transmission of pathogens with the administration of the anti-D prophylaxis (human immunoglobulin)
- Allergoid/anaphylactic reactions
- Shortened inpatient hospitalization for RhD-negative mothers since they can receive anti-D prophylaxis immediately after delivery if the test result is positive
- The anti-D immunoglobulin is currently obtained from immunized donors and imported



Limits of the test

- False-positive test results are possible in the case of a rare genotype, non-functional RhD variant¹ or a bone marrow donation⁴
- A test result cannot be reported if the mother is a carrier of the RhD gene
- False-negative results are also fundamentally possible^{5,6}

